



Liquid

Analysis



Systems

Components



Services

# Technical Information

# Oxymax H COS21

Sensor with long-term stability for frequent sterilization and autoclaving



#### Application

- Process control in enzyme production
- Control of culture growth
- Biotechnological production
- Food industry
- General process applications

#### Your benefits

- Sensor version suitable for the pharmaceutical industry:
  - Stainless steel 1.4435 (AISI 316L)
- Sterilizable and autoclavable
- Application-specific versions:
   Sensor for standard applications
- Sensor for standard
   Versatile usage:
  - Standard process connection Pg 13.5
  - Installation in standard pH assemblies possible
- Short response time:  $t_{98} < 60$  s
- Integrated temperature sensor



### Function and system design

Measuring	principle

The oxygen molecules diffused through the membrane are reduced to hydroxide ions (OH-) at the cathode. Silver is oxidized to silver ions (Ag+) at the anode (this forms a silver halogenide layer). A current flows due to the electron donation at the cathode and the electron acceptance at the anode. Under constant conditions, this flow is proportional to the oxygen content of the medium. This current is converted in the transmitter and indicated on the display as an oxygen concentration in mg/l, as a saturation index in % SAT or as an oxygen partial pressure in hPa.

#### Measuring system

- A complete measuring system comprises:
- Oxygen sensor COS21
- Transmitter, e.g. Liquisys M COM223/253 F
- Special measuring cable COK21
- Assembly, e.g. retractable assembly CPA475



- Fig. 1: Measuring system (example)
- 1 Oxygen sensor COS21
- 2 Retractable assembly CPA475
- *3* Special measuring cable COK21
- 4 Transmitter Liquisys M COM253 F

#### Input

Measured variable	dissolved oxygen [mg/l / % SAT / hPa]	
Measuring range	0,01 to 20 mg/l (ppm) 0 to 200% SAT	

### Wiring

**Electrical connection** 

The electrical connection between the sensor and the transmitter uses the multi-core special measuring cable COK21.



### Performance characteristics

Response time	From air to nitrogen at 25 °C (77 °F) • $t_{90} : < 30 s$ • $t_{98} : < 60 s$		
Reference operating conditions	Reference temperature: Reference pressure:	25 °C (77 °F) 1013 hPa (15 psi)	
Minimum flow rate	0.03 m/s (0.1 ft/s)		
Zero current	< 0.1 % of the current at air		
Maximum measured error	$\pm 1$ % of measured value <sup>1)</sup>		
Repeatability	$\pm 0.1$ % of measuring range end		
Long-term drift	Zero-point drift: Measuring range drift:	<0.1 % per week at 30 °C (86 °F) and under constant conditions $<0.1$ % per week at 30 °C (86 °F) and under constant conditions	
Influence of medium temperature	approx. 3.1% per K (compensated by transmitter)		
Influence of medium pressure	Pressure compensation not necessary		
Polarization time	< 60 minutes		
Polarization voltage	$-670 \pm 50 \text{ mV}$		
Oxygen intrinsic consumption	approx. 20 ng/h in air at 25 °C (77 °F)		

<sup>1)</sup> In accordance with IEC 746-1 at nominal operating conditions

# Installation

Angle of installation



Permitted installation angle

# Environment

Ambient temperature range	0 to +130 °C (32 to 270 °F)
Storage temperature	-10 to +60 °C (10 to 140 °F) at 95% relative air humidity, not condensing Caution! Danger of drying out Only store the sensor with the electrode protection cap (filled with 0.02 n NaOH).

		Process
Process temperature		0 to 130 °C (32 to 270 °F) Note! Please note the temperature compensation range of the transmitter used.
Process pressure		0 to 4 bar (0 to 58 psi)
Sterilization		When autoclaving, tightly place protection cap on plug-in head. No plug-in head protection is required when the device is sterilised in-situ.
		Note! If water has penetrated the plug-in head, dry it with warm or compressed air. This will prevent corrosion or contact problems.

## Mechanical construction

Design, dimensions	
	I       Threaded connection Pg 13.5         2,3       Seals         4       Cathode         5       Protective cap         6       Membrane         7       Membrane         7       Membrane         6       Membrane         7       Membrane         7
Weight	0.2 kg (0.44 lbs)
Material	Sensor shaft:Stainless steel 1.4435 (AISI 316L)Membrane:Optiflow®Electrode combination:Silver / PlatinumSealing ring:Viton®
Process connection	Thread Pg 13.5
Temperature sensor	NTC 22 k $\Omega$
Electrolyte	Alkaline electrolyte

# Ordering information

Product structure	I	Immersion depth	
	1	120 mm (4.72 in)	
	2	225 mm (8.86 in)	
	3	350 mm (13.8 in)	
	5	120 mm (4.72 in) + EN10204-3.1	
	6	225 mm (8.86 in) + EN10204-3.1	
	7	350 mm (13.8 in) + EN10204-3.1	
		Head	
		K Process Pg 13.5, four-pole plug-in head	
		Application	
		0 Hygienic version	
	COS21-	Complete order code	

Scope of delivery	<ul> <li>The following items are included in the delivery:</li> <li>Oxygen sensor with transport protection cap for membrane protection</li> <li>Electrolyte, 1 bottle, 50 ml (1.7 fl.oz.)</li> <li>Pipette for filling with electrolyte</li> </ul>		
Ì	Accessories Note! In the following sections, you find the accessories available at the time of issue of this documentation. For information on accessories that are not listed here, please contact your responsible service.		
Assemblies	Cleanfit H CPA475 • Retractable assembly for installation in tanks and pipework under sterile conditions • Technical Information TI240/C/07/en Unifit H CPA442		
	<ul> <li>Installation assembly for food, biotechnology and pharmaceuticals, with EHEDG and 3A certificate</li> <li>Technical Information TI306/C/07/en</li> </ul>		
Zero solution	<ul><li>3 units to produce 3 x 1 liter oxygen-free solution</li><li>order no. 50001041</li></ul>		
Electrolyte solutions and membrane cap kits	Electrolyte solutions ■ order no. 51505873		
	<ul> <li>Membrane kits</li> <li>Membrane kit Standard, COS21/COS21D: <ul> <li>order no. 51505874</li> </ul> </li> <li>Membrane kit Standard, COS21/COS21D, EN10204: <ul> <li>order no. 51516339</li> </ul> </li> <li>Membrane kit CIP, COS21/COS21D: <ul> <li>order no. 51518699</li> </ul> </li> <li>Membrane kit CIP, COS21/COS21D, EN10204: <ul> <li>order no. 71023225</li> </ul> </li> <li>Membrane kit FDA, COS21/COS21D: <ul> <li>order no. 71003199</li> </ul> </li> <li>Membrane kit FDA, COS21/COS21D, EN10204: <ul> <li>order no. 71023226</li> </ul> </li> </ul>		
	<ul> <li>Scope of delivery (all kits):</li> <li>3 Membrane caps</li> <li>1 O-ring (process seal)</li> <li>1 O-ring (sensor)</li> </ul>		
Measuring cable	<ul> <li>Special measuring cable for COS21</li> <li>COK21; cable length 3 m (9.8 ft) order no. 51505870</li> <li>COK21; cable length 10 m (33 ft) order no. 51505868</li> </ul>		
Transmitter	<ul> <li>Liquisys M COM 223/253 F</li> <li>Transmitter for oxygen, field or panel-mounted housing</li> <li>Hart<sup>®</sup> or PROFIBUS available,</li> <li>Ordering acc. to product structure, see Technical Information</li> </ul>		

### Documentation

- Technical Information Liquisys M COM223/253 F, TI246C/07/en
  Technical Information Cleanfit CPA475, TI240C/07/en
  Technical Information Unifit CPA442, TI297C/07/en

#### International Headquarters

Endress+Hauser GmbH+Co. KG Instruments International Colmarer Str. 6 79576 Weil am Rhein Deutschland

Tel. +49 76 21 9 75 02 Fax +49 76 21 9 75 34 5 www.endress.com info@ii.endress.com

TI244C/07/en/09.06 51505857 Printed in Germany / FM+SGML 6.0 / DT



People for Process Automation